ACNE LIGHT

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INDICATION AND FUNCTION OF EQUIPMENT

Acne Light has the function of emitting blue light for the treatment of acne in the following stages:

- Acne in the initial stage, slightly reddish;
- Papulopustulous acne with little inflammation.

The following situations are not handled by the Acne Light:

- Blackheads and milia not reddish and not inflamed;
- Papulopustulous acne followed by severe inflammation.

CLASSIFICATION

- **Standard:** RDC 185/2001 (ANVISA)
- **Classification:** II

SPECIFICATIONS

General Characteristics

Characteristics	Specifications
Power Supply	Supply Voltage: 90-240V~
	Current type: Alternate current
	Input Power: 5 VA
	Supply Frequency: 50/60Hz
	Operating Mode: Continuous operation
Protection to Water Penetration	IPX0
Dimension	27 x 88 mm (diameter x height)
Weight	60 g (± 1)
Treatment Time	2 minutes and 30 seconds

Autonomy Time of Battery	1 hour and 40 minutes
Battery	Lithium Ion Battery
	Maximum time of 1 hour to load battery

SAFETY - IMPORTANT PRECAUTIONS


-  Never look directly into the emitted LED light and especially do not direct it to any person unless that under treatment.
-  Only the components mentioned herein may be used in conjunction with the equipment.
-  In case any side effect is observed due to the use of the product, stop treatment and consult an esthetic professional.


-  Keep this product away from heat and water. The exposure to heat or water may cause malfunction.
-  Do not use the product around the eyes, lips, or directly on the eyes.
-  This product should not be used by children without adult supervision.
-  Only professionals of the esthetic sector should use this product (doctors, biomedicians, physiotherapists, technical and postgraduate estheticians).
-  Do not use the equipment concomitantly with isotretinoin-based drugs.

LIST OF COMPONENTS


- 01 Acne Light;
- 01 Power Supply;
- 01 User Guide;
- 01 Warranty Term.

PARAMETERS OF EQUIPMENT



OPERATION INSTRUCTION

To start treatment, quickly press the on/off button (2). Two quick sound signals will be emitted, the light indicator (1) will turn to green and the blue light will be emitted.

During treatment, a sound signal will be emitted every 10 seconds. Press the button (2) again at any time to stop or continue the treatment.

During treatment, the light beam should be directed to the desired area (slightly touching). At the end of the application, two sound signals will be emitted and the equipment will automatically turn off.

To turn off before concluding the treatment, keep the button pressed until two short sound signals are emitted after a long sound signal.

The skin may become temporarily red on the treated region.

The cover should be removed during application and placed back whenever the equipment is not being used.

The emission of four sound signals during treatment indicates that the battery of the equipment is weak and must be turned off.

It is recommended to use the product 2-3 times a day with an interval of three hours between applications.

After treatment, the skin may become red, which normally disappears within 1 hour. If the skin does not recover after a long time, interrupt the use of the product and seek an esthetician.

Only use the product after charging it for at least 30 minutes. The tip for the treatment must be disinfected before use.

CLEANING OF EQUIPMENT

To hygienize the product, the treatment pointer must be cleaned with a clean cloth moistened with 70% alcohol.

The product is NOT AUTOCLAVABLE. If it is exposed to that process, the LED emitter will suffer permanent damage.

CHARGING ACNE LIGHT

- Connect the plug to the power supply on the terminal located on the lower part of the Acne Light;
- The light indicator may present the following:
 - Blinks green and red: Battery charging;
 - Continuous green: Charged battery.
- With the product completely charged, it may be carried out over 40 treatments.

For safety reasons, do not operate the product while it is connected to the charger.

PREVENTIVE MAINTENANCE

The equipment must be calibrated at least every two years by the manufacturer. If this maintenance is not performed, the manufacturer assumes no responsibility for the safe operation of equipment.

All support services, such as alterations, repairs, calibration, etc., can only be performed by the manufacturer. The circuit schematics, parts lists, descriptions, instructions for calibration and measurement are not available to persons not qualified by the manufacturer.

If maintenance or other support service is performed by unauthorized service personnel, the manufacturer assumes no responsibility for the safe operation of equipment.

Never open the equipment. Any problems call the technical support of DMC Equipamentos.

STORAGE AND TRANSPORTATION

- Store the equipment protected from dust, exposure to direct sunlight, away from chemicals and cleaning agents;
- The equipment should be stored, transported and used in the following environmental conditions:
 - Temperature: +10 °C to +40 °C;
 - Humidity: 30 % to 75 %;
 - Atmosphere Pressure: 700 hPa to 1060 hPa.

Avoid equipment fall.

Keep the equipment in a safe place, avoiding subjecting him to blows and vibration.

DISPOSAL

After the lifetime of the product and its accessories, it can cause environmental contamination or may be misused. To minimize these risks, the client should discard the product as determined by local law.

TROUBLESHOOTING

Type of error	Possible Solutions
The equipment does not turn on.	Charge the battery. If after these steps the machine still has problems, send to the manufacturer's technical support.
Does not emit blue light	Send the equipment to the manufacturer's technical support.

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Proteja contra radiação
*Proteja contra radiación
Protect it from radiation*

Manter seco
*Mantener seco
Keep it dry*

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A. The equipment manufactured and / or marketed by the DMC are guaranteed for 24 (twenty four) months from the date of customer purchase, against manufacturing defects.

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B. The warranty only covers defects in manufacturing or materials used to in manufacturing of products. Warranty does not cover shipping expenses.

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C. The warranty is automatically canceled if occur electrical, physical abuses, if the parts are changed, or if there are applications other than those for which the equipment was developed.

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D. For equipment repaired out of warranty period, it will only be extended to components replaced.

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E. The most common causes of defects are due to physical shock applied to the device, where the warranty is canceled.

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F. The DMC is not responsible for damage or injury resulting from misuse of equipment it produced and / or sold, leaving it to the user provide security measures to prevent such occurrences.

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G. The responsibility of the DMC regarding the use of the equipment and its consequences is limited to the replacement value thereof.

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Congratulations! You have just purchased a high quality DMC product. The next step is to validate the warranty of your equipment. DMC provides you with two options:

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Access the electronic address of DMC www.dmcgroup.com.br click on PRODUCTS > REGISTRATION OF PRODUCTS. Fill out the form and conclude the registration. And it is ready!

The warranty of your product will already be validated.

CORREIO / REPRESENTANTE REPRESENTATIVE

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Fill out this form (page 34), detach it and send to a DMC Representative of your country.

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